

# GyroWash

---

Available as 8 and 20 pot instruments, the GyroWash is a colour fastness to laundering tester for textiles and leather, according to AATCC and ISO standards.

Featuring the most simple to use touchscreen user interface on the market.



---



# GyroWash

## COLOUR FASTNESS TO LAUNDERING TESTER

GyroWash colour fastness to laundering tester features the latest state of the art touchscreen. Designed alongside laboratory technicians, the new touchscreen user interface features pre-loaded laundering standards to aid the technician with setting up the test.



### KEY BENEFITS

#### TOUCHSCREEN USER INTERFACE

An intuitive touchscreen user interface, designed alongside textile technicians in a working laboratory, offers testing which is simple to set up. Pre-loaded standards, temperature settings, rotational speed and an autostart function can all be controlled from the touchscreen.

#### CHOICE OF CAPACITY

The GyroWash comes in both 8 and 20 pot models, with a choice of 1200ml and 500ml pots.

#### INTERCHANGEABLE VESSELS

A push and rotate twist motion makes changing between 1200ml and 500ml pots quick and easy.

#### COOL TOUCH LID

Safe access to samples after the instrument has operated at extremely high temperatures is made possible through a cool touch lid.

#### INSULATED BATH

A sealed and insulated bath results in lower energy usage, making the instrument more efficient.

#### HIGH WASHING TEMPERATURE LIDS

For textile testing laboratories performing colour fastness to washing and cleansing procedures above 60°C, James Heal have developed high washing temperature lids which contain the wash liquor without leakage. These guarantee complete confidence in results, and accessories allow for easy access after the test.

#### A RANGE OF STANDARDS

Both ISO and AATCC standards, and several retailer methods, can be tested on GyroWash.

#### QUICK DRAIN AND FILL

Easy to reach fittings on the front of the instrument make draining and filling the instrument with water easy to do.

#### DEDICATED PREPARATION AREA

An integrated drip tray collects any surplus water which drains away without the need for emptying.

#### FULL RANGE OF TEST MATERIALS

A comprehensive selection of high quality test materials, including greyscales, detergents and multifibres, are available for the GyroWash.

Dimensions (mm)	Height (mm)	Width (mm)	Depth (mm)	Weight (kg)
GyroWash - 8 pot	1040	876	731	135.0
GyroWash - 20 pot	1040	1266	731	238.0

### CONTACT US

#### For more information

web: [www.james-heal.co.uk/en/gyrowash](http://www.james-heal.co.uk/en/gyrowash)

tel: +44 (0) 1422 366 355

email: [sales@james-heal.co.uk](mailto:sales@james-heal.co.uk)

James Heal Richmond Works Halifax UK HX3 6EP

 **James Heal**  
Extraordinary Testing Solutions